



## WRITING STAND WITH 1 DRAWER AND 2 SHELVES ON CASTORS



<b>Length</b>	550
<b>Width</b>	500, 600, 700

Frame structure in stainless steel AISI 304 fully welded 40x40mm tube, 1.00mm thick. With one drawer mounted on stainless steel rails. With 2 shelves manufactured in stainless steel AISI 304, welded to the table framework. Inclined writing top with front bended edge to prevent paper drop. Frontal handle in round pipe. On 2 swivels and 2 braked anti-trace castors, diam. 125mm, and with stainless steel frame.

---

**SKU:** N/A



		W		
		500	600	700
L	Cod.	Cod.	Cod.	
	550	T10501001	T10502001	T10503001

**W** = Width / Profondità / Breite / Profondeur.

**L** = Length / Larghezza / Länge / Longueur.