



WRITING STAND WITH 1 DRAWER AND 2 SHELVES ON CASTORS



Length	550
Width	500, 600, 700

Frame structure in stainless steel AISI 304 fully welded 40x40mm tube, 1.00mm thick. With one drawer mounted on stainless steel rails. With 2 shelves manufactured in stainless steel AISI 304, welded to the table framework. Inclined writing top with front bended edge to prevent paper drop. Frontal handle in round pipe. On 2 swivels and 2 braked anti-trace castors, diam. 125mm, and with stainless steel frame.

SKU: N/A



		W		
		500	600	700
L	Cod.	Cod.	Cod.	
	550	T10501001	T10502001	T10503001

W = Width / Profondità / Breite / Profondeur.

L = Length / Larghezza / Länge / Longueur.