



MOBILE UNITS WITH 4 DRAWERS




Length	500, 700
Width	600, 700, 800

Manufactured in stainless steel AISI 304, with scotch brite dull finishing. Full extraction drawers mounted on stainless steel rails. Upper top with sandwich structure, with front grip in stainless steel round wire. On 2 swivels and 2 braked anti-trace castors, diam. 125mm, and with stainless steel frame.

SKU: N/A



L	W		
	600	700	800
Cod.	Cod.	Cod.	Cod.



500	T04201001	T04202001	T04203001
700	T04201002	T04202002	T04203002

W = Width / Profondità / Breite / Profondeur.
L = Length / Larghezza / Länge / Longueur.