



PERFORATED TABLES FOR LAMINAR FLOW WITH LOWER SHELF



Length	500, 600, 700, 800, 900, 1000, 1100, 1200, 1300, 1400, 1500, 1600, 1700, 1800, 1900, 2000
Depth	600, 700, 800

Frame structure in stainless steel AISI 304 fully welded 40x40mm tube, 1.00mm thick, rear legs 80mm advanced, suitable for rounded floor corner. Rear and side welded reinforcement rails. Perforated working top in AISI 304 sheet, 1.50mm thick. Hole diameter 9mm and perimetral frame 60mm wide. Rounded (R=10mm) front upper edge, as well as rounded (R=2mm) side and back edges. Lower perforated shelf fully welded to the framework. On adjustable feet.

SKU: N/A



L	W		
	600	700	800
	Cod.	Cod.	Cod.
500	T06101001	T06102001	T06103001
600	T06101002	T06102002	T06103002
700	T06101003	T06102003	T06103003
800	T06101004	T06102004	T06103004
900	T06101005	T06102005	T06103005
1000	T06101006	T06102006	T06103006
1100	T06101007	T06102007	T06103007
1200	T06101008	T06102008	T06103008
1300	T06101009	T06102009	T06103009
1400	T06101010	T06102010	T06103010
1500	T06101011	T06102011	T06103011
1600	T06101012	T06102012	T06103012
1700	T06101013	T06102013	T06103013
1800	T06101014	T06102014	T06103014
1900	T06101015	T06102015	T06103015
2000	T06101016	T06102016	T06103016

W = Width / Profondità / Breite / Profondeur.

L = Length / Larghezza / Länge / Longueur.