



SERVICE TROLLEY WITH 2 SCREWED PERFORATED SHELVES ON ANTI-TRACE CASTORS



Length	500, 600
Width	600, 800, 900, 1000

Trolley manufactured in stainless steel AISI 304. Frame in round tube with 2 screwed perforated shelves H=40mm, equipped with perimetric containment rim. With pull handle. On 2 swivels and 2 braked anti-trace castors, diam. 125mm, and with stainless steel frame.

SKU: N/A



		W			
		600	800	900	1000
L	Cod.	Cod.	Cod.	Cod.	Cod.
500	T93601001	T93602001	T93603001	T93604001	
600	-	-	T93603002	T93604002	

W = Width / Profondità / Breite / Profondeur.

L = Length / Larghezza / Länge / Longueur.