



## WEIGHING TABLE WITH CENTRAL GRANITE PLATE



<b>Length</b>	500, 600, 700
<b>Width</b>	600, 700, 800


Frame structure in stainless steel AISI 304 fully welded 40x40mm tube, 1.00mm thick, rear legs 80mm advanced, suitable for rounded floor corner. Rear and side welded reinforcement rails. Working top in AISI 304 sheet, 1.50mm thick, with inner plywood reinforcement and stainless steel lower cover. Opening/s on the work surface to allow the housing of 1 or 2 400x400x40Hmm black granite plates, supported by anti-vibration rubber spacers. Rounded (R=10mm) front upper edge, as well as rounded (R=2mm) side and back edges. On adjustable feet.

---

**SKU:** N/A



		W		
		600	700	800
L	Cod.	Cod.	Cod.	Cod.



500	T07001001	T07002001	T07003001
600	T07001002	T07002002	T07003002
700	T07001003	T07002003	T07003003

**W** = Width / Profondità / Breite / Profondeur.

**L** = Length / Larghezza / Länge / Longueur.