




WEIGHING TABLE WITH DOUBLE GRANITE PLATES

Length 1400, 1600, 1800

Width 600, 700, 800



Frame structure in stainless steel AISI 304 fully welded 40x40mm tube, 1.00mm thick, rear legs 80mm advanced, suitable for rounded floor corner. Rear and side welded reinforcement rails. Working top in AISI 304 sheet, 1.50mm thick, with inner plywood reinforcement and stainless steel lower cover. Opening/s on the work surface to allow the housing of 1 or 2 400x400x40Hmm black granite plates, supported by anti-vibration rubber spacers. Rounded (R=10mm) front upper edge, as well as rounded (R=2mm) side and back edges. On adjustable feet.

SKU: N/A



		W		
		600	700	800
L	Cod.	Cod.	Cod.	Cod.

	1400	T07301001	T07302001	T07303001
	1600	T07301002	T07302002	T07303002
	1800	T07301003	T07302003	T07303003

W = Width / Profondità / Breite / Profondeur.

L = Length / Larghezza / Länge / Longueur.